



Instructions

What's included:

- Vi-band head massager
- Micro USB charging cable

Operation:

- Charge Vi-Band fully (see charging)
- Align **center message node** over middle of forehead
- Align **outer message nodes** over your temples
- Strap velcro so that the band is **snug** against your forehead
- Power **ON** device using **slide switch** on the right side of the case.
- **Blue LED** illuminates when power is **ON**.
- Cycle massage pattern by pressing **mode button** in center of the case.

(Flip over)

Instructions Continued

Charging:

1. Insert **USB connector** of **charging cable** into any compatible USB charger.
2. Plug **micro USB connector** of **charging cable** into **micro USB port** on left side of case. **Charge LED** will illuminate.
3. After about 2 hours, the **charge LED** will turn off, indicating that the unit is **fully charged**.

Cleaning:

- The Vi-Band is **surface clean** or **dry clean only**
- The Vi-Band is **not waterproof**; please do not submerge in water

For more information, please visit vi-band.com

