

What's so wonderful about
Intellinetix® Vibrating Calf/Shin Wrap?



Easy to Use • Easy to Care For • Easy to Maintain • Easy to Believe In

Easy to **maintain.**

To recharge, simply plug one end of charger into an AC wall outlet and the other end into the charge port on the Vibration Pad.

A green light will illuminate when the unit is done charging.

Product box is perfect for storage or travel.

Easy to **service.**

Brownmed warrants this product to be free from manufacturing defects.

If the product stops working within one year from date of purchase, we will replace it.

Need a replacement or extra parts? Customer service is happy to help you. Please call us at **1-800-843-4395** or **+1-712-336-4395**.

Vibration Pad reorder number: 07516
Charger reorder number: RM07502

Easy to **care for.**

To clean, remove Vibration Pad and machine wash Calf/Shin Wrap in cold water. Lay flat and allow to air dry. Spot clean Vibration Pad as needed.

Keep charger opening dry while cleaning.

Caution:

Do not allow the Vibration Pad to get wet.

For adult use only.

Do not dry clean.



Easy to **use**.

- 1 Before first use, charge Vibrating Calf/Shin Wrap for four full hours. Fully charged, Calf/Shin Wrap will operate up to 30 minutes.
- 2 Slip the device on and pull up over your calf and shin.
- 3 Ensure Vibration Pad covers source of pain.
- 4 Apply vibration therapy for 25-30 minutes for maximal duration of pain relief.
Note: Vibrating Calf/Shin Wrap can be turned on and off at your leisure by pressing the button near the top of the Vibration Pad.
- 5 Replenish battery after each use by plugging charger into the charge port on the removable Vibration Pad.

Easy to **believe in.**

Intellinetix® Guiding Principles

Use **innovation** to simplify traditional treatments.

Have an **optimistic** attitude that begins with a desire to understand and respond to patient needs.

Believe **simple is best.**

Know **inspiration** can happen anywhere.

Encourage **wellness** through positive affirmations.

Elevate the user experience with **words of encouragement** and **healthy living suggestions** that are tied into each Intellinetix product and package.

Love the Vibrating Calf/Shin Wrap and want to know what other products are in the Intellinetix line?

Visit us on the web at **www.brownmed.com**



Executive:

280 Summer Street, Suite 400
Boston, MA 02210
P 857-317-3354
F 857-317-3355

Sales:

4435 Main Street, Suite 820
Kansas City, MO 64111
TF 877-853-5518
P 816-581-7001
F 816-581-7005

Manufacturing:

1300 Lundberg Drive West
Spirit Lake, IA 51360
TF 800-843-4395
P 712-336-4395
F 712-336-2874