

SHARPER IMAGE®

T.E.N.S. FOOT MASSAGER WITH INFRARED HEAT

Item No. 205748

User Guide



Thank you for purchasing the Sharper Image T.E.N.S. Foot Massager with Infrared Heat. Please read this guide and store it for future reference.

INTRODUCTION

The foot is one of the most sensitive parts of the body, with a complex network of muscles, nerves and blood vessels. We depend on our feet every day, so it is only natural that we would want to keep them relaxed, refreshed and restored. In fact, many practitioners of foot reflexology believe that the health of the foot is directly related to the health of the heart, the liver, the kidneys and other vital organs.

The T.E.N.S. Foot Massager with Infrared Heat is designed to provide temporary relief from stress, fatigue and soreness in the feet. As an added feature, you may use the included electrode cables, belt and gel pads to administer T.E.N.S. massage to other parts of the body, including the legs, arms, back, etc. In addition to T.E.N.S. massage, this device provides Far Infrared Ray [FIR] heat for additional relief.

WHAT IS T.E.N.S.?

Transcutaneous Electrical Nerve Stimulation [T.E.N.S.] is the advanced electrotherapy used by doctors and chiropractors worldwide. T.E.N.S. is a therapy that uses tiny electrical impulses that block pain signals before they can reach the brain. These gentle impulses also promote blood flow, which loosens and relaxes tight muscles. Finally, T.E.N.S. stimulates the production of endorphins, the body's own natural pain killers.


Several years ago, Sharper Image revolutionized the massage industry by adapting professional T.E.N.S. therapy for the home market in a convenient, portable size. The Sharper Image T.E.N.S. Foot Massager with Infrared Heat has been designed and produced to our high standards, in order to deliver soothing and invigorating relief for feet and other parts of the body.

FEATURES AND BENEFITS:

- T.E.N.S. and FIR heat can be used together or independently
- 25 Pulse Modes
- Automatic shutoff timer [10-90 minutes]
- Includes Electrode cables and self-adhesive gel pads provide T.E.N.S. massage for other parts of the body
- Includes electrode belt for T.E.N.S. massage on the lower back or abdomen
- Low frequency electrical pulses stimulate the muscles, promote blood flow and the production of endorphins [the body's own natural pain killers]
- Far Infrared Ray [FIR] heat penetrates and soothes muscle and bone tissue without feeling hot to the touch
- Drug-free, non-invasive relief for stress, fatigue and sore muscles
- Multiple levels of electromagnetic wave intensity

UNIT COMPONENTS

1. Far Infrared Ray [FIR] heat button
2. Pulse electromagnetic wave display screen
3. Timer display screen
4. Mode switch display screen
5. Body intensity display screen
6. Soles intensity display screen
7. Electrode pads for soles
8. On/Off button
9. Increase value function button [+]
10. Decrease value function button [-]
11. Function button
12. Auto/Mode switch button
13. Electrode jack
14. DC power jack
15. Infrared heat zones



INCLUDED ACCESSORIES

1. Electrode cables
2. Self-adhesive electrode gel pads
3. Electrode belt
4. AC power adapter
5. Wireless remote control



USING YOUR T.E.N.S. FOOT MASSAGER WITH INFRARED HEAT:

1. Place unit on a flat, stable surface.
2. Plug unit into an AC outlet.
3. Press the POWER button to turn unit on.
NOTE: When the device is powered on, the default setting is: TIME = 25 minutes, AUTO/MODE = automatic, SOLE = 0 intensity/mode and BODY = 0 intensity/mode. The SOLE display will be flashing. Press “+” or “-” to set the massage mode and intensity.
4. Use the FUNCTION button to cycle through the four different functions: TIME, AUTO/MODE, BODY and SOLE.
5. Place your clean, dry, bare feet on the foot pads.
6. Select SOLE and press “+” or “-” to increase or decrease the intensity/mode to the desired level.
7. Select MODE and press “+” or “-” to select the desired MODE setting.

PLEASE NOTE: BOTH feet must be bare and on unit to feel electrical pulses.

USING YOUR T.E.N.S. MASSAGER FUNCTION WITH YOUR SELF-ADHESIVE ELECTRODE GEL PADS AND ELECTRODE BELT:

1. Connect the electrode cables to the self-adhesive gel pads or the electrode belt.
2. Apply the self-adhesive gel pads to clean skin, or wrap the electrode belt around your waist.
3. Plug the electrode cables into the electrode jack.
4. Select BODY and press “+” or “-” to increase or decrease the intensity to the desired level.
5. Select MODE and press “+” or “-” to select the desired MODE setting.

INTENSITY/PULSE MODES

Mode #	Description
1	Intermittent soft beat
2	Continuous powerful beat
3	Right and left continuous powerful beat
4	Asynchronous gentle tap
5	Right and left heavy massage
6	Light tap
7	Hard tap
8	Massage
9	Hammer
10	Beat
11	Light beat
12	Invigorating wave

13	Right and left jog
14	Light press
15	Hard hustle
16	Hard tingle
17	Vibration
18	Intense knead
19	Right and left heavy pulse
20	Right and left heavy beat
21	Right and left gentle hammer
22	Rhythmic tap
23	Powerful pulse
24	Powerful hammer
25	Right and left heavy press

SAFETY WARNINGS

- Do not use this device near water. Do not immerse in water or any other liquid.
- Keep out of reach of children at all times.
- Clean this unit with a soft, dry cloth when needed.
- Always start with a low intensity level, then and work your way up to a stronger intensity.
- Stop using this device if you feel discomfort or nausea, or experience redness.
- Not for use if you have a pacemaker, defibrillator or other metallic device.
- Not for use if you have epilepsy, cancerous lesions, tumors, acute trauma, fractures, hernias, critical ischemia of lower limbs, blood flow deficiencies, venous thrombosis, or if you are recovering from surgery.
- Not for use by pregnant women.
- Unplug this device when not in use.
- Consult your doctor before using any at-home massage product.

WARRANTY/CUSTOMER SERVICE

Sharper Image branded items purchased from SharperImage.com include a 1-year limited replacement warranty. If you have any questions not covered in this guide, please call our Customer Service department at 1 (877) 210-3449. Customer Service agents are available Monday through Friday, 9:00 a.m. to 6:00 p.m. ET.

SHARPER IMAGE®

The SHARPER IMAGE® name and logo are registered trademarks.

©Sharper Image. All Rights Reserved.