



Manufactured by

LED Technologies, Inc.  
Greenwood Village, CO 80111  
(800) 337-9565 | Made In China  
revivelighttherapy.com

## CUSTOMER SERVICE

You can reach Customer Service by calling 1(800)337-9565, Monday-Friday, 8:00AM-4:30PM, Mountain Standard Time. Customer Service can also be reached by email at support@revivelighttherapy.com.

## REPLACEMENT PARTS

The following list of damaged parts can be purchased as replacement items directly from LED Technologies, Inc. by calling Customer Service at 1(800)337-9565:

- Power Adapter
- Treatment Brushes

## STAY CONNECTED

When you register your Soniqué™ system at revivelighttherapy.com/warranty, you will be entered into our system to be eligible for exclusive offers, VIP discounts, new product previews, and the opportunity to earn free products through new product trials.

## LIMITED WARRANTY INFORMATION

Congratulations on your purchase. We value you as a client and strive for the highest quality products to earn your satisfaction. Please use this as your reference for any customer service or warranty needs you may have. You must register your product within 30 days of purchase at www.revivelighttherapy.com for the warranty to be effective.

## WHAT IS COVERED BY THE LIMITED WARRANTY?

This limited warranty is effective only if the device is purchased and used in the United States. This 1 year limited warranty covers defects in materials and workmanship in your Soniqué™ LED sonic cleansing system, if you purchased it from an authorized partner or directly from the manufacturer (LED Technologies, Inc.). Your 1 year limited warranty is valid for you alone, the original purchaser, and cannot be transferred. If the product was a gift, you can ask the gift giver to send a copy of the receipt to support@revivelighttherapy.com, or process the warranty claim on your behalf.

## WHAT ISN'T COVERED BY THE LIMITED WARRANTY?

This limited warranty does not cover damage, defects, or failure caused by abuse or misuse, improper care, impact or insertion of foreign objects, use with incorrect power adapter, or loss of device.

## HOW LONG DOES THE LIMITED WARRANTY LAST?

The limited warranty begins from the date of purchase and expires 1 year later. The warranty period is not extended if your device is replaced by LED Technologies, Inc. during the warranty period.

## WHAT WILL WE DO IF YOUR DEVICE IS DEFECTIVE?

If you discover a covered defect in your Soniqué™ LED sonic cleansing system before the 1 year warranty period expires, contact Customer Service at 1(800)337-9565 for assistance. All limited warranty services must be authorized and performed by LED Technologies, Inc.

## WARRANTY LIMITATIONS & DISCLAIMERS

This warranty gives you specific legal rights, and is in addition to your statutory rights. Your remedies for the breach of this limited warranty are limited solely to the repair and replacement of your device as set forth in this warranty statement. To the maximum extent permitted by applicable law, LED Technologies, Inc. disclaims all implied warranties. If applicable law does not permit the disclaimer of implied warranties, then the duration of implied warranties are hereby expressly limited to the duration of the express limited warranty set forth above.

To the maximum extent permitted by applicable law, in no event shall LED Technologies, Inc. be liable for any special, incidental, consequential, punitive, or indirect damages based upon breach of warranty, breach of contract, negligence, strict product liability, tort or other legal theory, even if LED Technologies, Inc. is advised or aware that such damages might arise.

Such damages include, but are not limited to, loss of device, cost of any substitute equipment, the claims of third parties, or damage to property. If the applicable law does not allow limits on warranties, or on remedies for breach in certain transactions, the limits in this paragraph and the disclaimer or exceptions from warranty may not apply. This limitation does not apply to claims for death or personal injury caused by a defect in the device or any other liability that cannot be excluded or limited by law. These warranty terms, conditions, and limitations allocate the risks of product failure between LED Technologies, Inc. and you, the purchaser of Soniqué™ LED sonic cleansing system. This allocation is recognized by both parties and is reflected in the price of Soniqué™ LED sonic cleansing system. By using Soniqué™ LED sonic cleansing system you acknowledge that you have read these terms, conditions, and limitations, understand them and are bound by them.



*Soniqué*™

Wrinkles  
LED Sonic Cleansing System  
by reVive Light Therapy®


Soniqué™ LED Sonic Cleansing System- sonic cleansing device with LED light therapy for treatment of wrinkles.

## SYSTEM INCLUDES

- Soniqué™ LED Sonic Cleansing System
- 2 Interchangeable Silicone Treatment Heads
- USB Charging Cord
- Power Adapter
- User Guide



Cleansing Brush Head



Massage Brush Head

## PRODUCT REPRESENTATION

LED Technologies, Inc. is constantly improving its products with the latest technology. As a result, the actual product may differ from the product shown in images on the packaging and in this user guide.

## CONTRAINDICATIONS

- Do not use this device if your skin is light sensitive or you are currently using medication that may cause skin sensitivity to light. Doing so may result in photosensitive reaction.
- If unsure whether your medication causes sensitivity to light, confirm with your physician prior to using this device.
- Discontinue use immediately if you experienced any sensitivity to using this device.

## WARNINGS

- Do not use the device if power adapter has fallen into water or is damaged in any way. Doing so could cause electric shock or result in injury. Please contact Customer Service at 1(800)337-9565 if power adapter has been damaged in any way.
- Stop using device immediately if you think it is defective. Using a damaged device could result in injury. Contact Customer Service at 1(800)337-9565.

- Do not bend or abuse the wire of the power adapter as this can cause electrical shock and problems.
- Use only the supplied power adapter. Use of incorrect power adapter could damage the device. Contact Customer Service at 1(800)337-9565.
- Do not attempt to take the device or power adapter apart. Any modification or attempted service will void the warranty.
- Do not store device in sun or location with high temperatures as this may cause damage to plastic casing.
- Discontinue use and contact customer service if device overheats or becomes uncomfortably hot to your skin.
- Do not exceed recommended treatment times or usage.

## PRECAUTIONS

- Always unplug power adapter from wall outlet when not in use.
- Do not use this device on infants or young children without physician approval.
- Do not shine the device lights directly into the eyes.
- This product contains medical grade silicone. If you have or suspect you have allergies to silicone, consult your Doctor prior to use.
- Contact your Doctor if you have any medical questions regarding use of this device.
- Contact your Doctor prior to using if you:
  - Are pregnant or suspect you may be pregnant.
  - Have suspicious or cancerous lesions.

## DISCLAIMER

This product is not intended to diagnose, cure or prevent any disease. LED Technologies, Inc. makes no claims, representations or warranties regarding the ability of this product to cure any physical, skin or mental condition. A qualified health professional should be consulted with regards to any condition requiring medical attention.

This product contains a lithium-ion battery that must be disposed of properly, according to local laws and regulations, and not as household waste. When the product reaches its end of life, check with your local state and municipal waste disposal regulations, to learn about disposal and recycling options.

## Soniqué™ LED Sonic Cleansing System

This is a skin cleansing device with light emitting diodes (LEDs) that emits energy for use in dermatology, for the treatment of fine lines, wrinkles, and other effects of the aging skin. The device delivers natural light energy in Infrared, Deep Red, Red, and Amber spectrum's that are proven to reduce and/or reverse fine lines and wrinkles, enhance the skin tone and texture, and create more youthful appearance. There are no user settings or adjustments.

## CHARGING INSTRUCTIONS

1. Unpack your device.
2. Plug USB charging cord into port connector on the device.
3. Plug USB charging cord into a USB power adapter, and into outlet for charging.
4. Charge your Soniqué™ LED sonic cleansing system for 2-4 hours. Device is fully charged when indicator light is solid.

## LIGHT SENSITIVITY TEST

Test yourself for light sensitivity prior to initial treatment:

1. Power on your device.
2. Place device on or within 1/4" of the skin on your forearm for three minutes.
3. After 3 minutes check the area. If area has turned red and lasts for more than 2 hours your skin is light sensitive.
4. Discontinue use and contact customer service.

## OPERATING INSTRUCTIONS FOR SONIC CLEANSING MODE

1. Fully charge device (2-4 hours).
2. If you use makeup, remove eye makeup by hand.
3. Attach silicone cleansing brush head.
4. Wet your skin and the brush head with warm water.
5. Liberally apply cleanser to your skin or the brush head.
6. Power up device (lights will turn on).
7. Press vibration button.
8. Gently move device over the treatment area in small circular motion.
9. Turn off device when finished. The device will auto shut-off after 3 minutes.

## OPERATING INSTRUCTIONS FOR LIGHT TREATMENT ONLY

1. Fully charge device (2-4 hours)
2. Clean and dry area prior to use.
3. If you will be using serums prior to light therapy treatment, attach silicone massage head.
4. If you are not going to be using serums prior to light therapy treatment, use the device without silicone treatment head.
5. Turn on device.
6. Place device directly on the skin over the desired treatment area for 3 minutes.
7. The device will auto shut off after completion of 3 minute treatment.
8. Move the device to the next treatment area and repeat step 5.
9. Use device daily for best results. Do not exceed 3 treatments per day per treatment area.

## MAINTENANCE & CLEANING

Maintenance:

1. This is a medical device. It is important to keep it clean.
2. The device has no user serviceable parts.
3. All settings are set at the factory.
4. Charge device only with supplied USB power adapter.
5. Power indicator light will begin to flash when battery levels are low, indicating it is appropriate to begin charging.

Cleaning:

1. Turn device off and unplug from any power source.
2. Remove brush head, rinse and dry.
3. Wipe device surface with damp cloth or alcohol wipes.

## TROUBLESHOOTING

1. System is not turning on:
  - Is the device fully charged? (Power indicator light will be on and solid.)
  - If device still does not turn on, call Customer Service at 1(800) 337-9565.
2. System is not charging:
  - Is USB charging cord firmly inserted into device? (Is power indicator light turned on?)
  - Try different USB power adapter to ensure it is working properly.
  - Try different wall outlet to determine if there is power to the outlet.
  - If device still does not charge, call Customer Service at 1(800) 337-9565.
3. System shuts off abruptly, prior to completion of 3-minute treatment:
  - System is designed to shut off automatically after 3 minutes.
  - Make sure device is fully charged (see charging instructions.)
  - If the device continues to shut off prior to 3 minutes, please contact Customer Service at 1(800) 337-9567
4. Some lights appear turned off:
  - The infrared (IR) light spectrum is not visible to the human eye. These lights will look as if they are not on but are operating properly. If the red visible lights are on, then the device is operating as designed.

## TRAVELING WITH DEVICE OUTSIDE USA

Power adapter is designed for USA standard 120V power outlets. When traveling outside the USA, you will need to check the local power requirements and use the appropriate power adapter and/or converter.